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* Easy Wiping Kit (BPW)

Gamma sterilization

* Buffered Peptone Water (BPW, ISO composition)*1 is prepared in accordance with the usage and is aseptically packaged.

Formula (in 20.07g for 1 liter)

Enzymatic digest of casein ······	10.00 g
Sodium chloride ······	5.00 g
Potassium dihydrogen phosphate ······	1.50 g
Disodium phosphate (anhydrous) ······	3.57 g
(Dodecahydrate of 9.00 g in anhydrous conversion)	
pH 7.0±0.2	


* *1 AccuDia™ Buffered Peptone Water (BPW, ISO composition) (Product code 05121), ISO 6887-1:2017 referenced

- The cotton swab cap of this product is orange.
- The product is sterilized by gamma irradiation after packaged in the case.

* ——— Directions ———


1) Pull out the orange-colored cotton swab cap. Using the cotton swab, wipe out a fixed amount of surface area of the inspection item.

Example) Surface area: For a flat surface, use the wipe frame 100 (code 06598) or equivalent to determine a fixed amount of surface area to wipe. In the case of an irregular, complex surface, verbalize the width, length and other measure of size in detail, and wipe the surface the same way each time.


Example) Wiping force to apply: Thrust the cotton swab onto the area of interest with a force that can almost bend its stick part (approximately 100 g to 150 g).

Example) How to wipe (100 cm²): Start wiping from the upper left corner laterally and go on to the bottom side, and then, start from the lower left corner vertically and go on to the right side, so that the rectangle can be wiped up evenly.


2) Put the cotton swab back in the body and seal it with the cap tightly. Tilt the body and mix well to make a specimen solution.

Example) Shake the wiping kit 20 times vigorously in an up-and-down reciprocating motion.

3) Remove the top cap on the edge, and push the center part of the body to inoculate all of the specimen solution (1 mL) onto Compact Dry™. As required, dilute the solution stepwise and inoculate the dilutions.

* ——— Precautions for use and handling ———

- 1) When pulling out the stick cap, be careful not to allow the cotton swab to touch your hand.
- 2) To maintain the wiping operation constant each time, we recommend that the operation method be verbalized by reference to the examples given in the operation steps in Directions.
- 3) Inspect the taken specimens within 30 minutes wherever possible. Where they are transported, condition them under between 1 and 8°C and have them inspected within 24 hours.
- 4) Easy Wiping Kit is coordinated to discharge 1 mL when squeezed out in the full amount.
- 5) For the inoculation, cultivation and judgment with Compact Dry™, follow the instruction manual of each Compact Dry™ product.

————— Storage and shelf life —————

[Storage]

Keep at room temperatures (1 ~ 30°C).

[Shelf life]

Two (2) years after the date of manufacture.

* ——— Package ———

Easy Wiping Kit (BPW) Gamma sterilization 1 mL × 200 units

····· Code 06537

* ——— Further information ———

Customer Support Section, Shimadzu Diagnostics Corporation
3-24-6, Ueno, Taito-ku, Tokyo 110-0005 Japan
TEL: +81-3-5846-5707

* Manufactured by
Shimadzu Diagnostics Corporation

3-24-6, Ueno, Taito-ku, Tokyo 110-0005 Japan