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(0186S)

**\*\* Test for *Escherichia coli* and coliforms in water**

**\*\* EC Blue-MPN plate**

The product is individually packed, and sterilized by  $\gamma$ -rays to be used directly after opening

**[Component]**

- EC Blue-MPN plate
- MPN plate main frame
- MPN plate lid
- MPN calibration scale

**[Directions]**

- ① EC Blue solution (100mL) is prepared using a watery sample and EC Blue 100/100P.
- ② The EC Blue solution is poured into the slope of an MPN plate. After each well (10, 1, 0.1mL) is filled, excess of the EC Blue solution is pooled into fluid reservoir.
- ③ The MPN plate is covered with a lid, then incubated at 35-37°C for 24 h ~28h.
- ④ After incubation, fluorescence is observed in wells by irradiating with a UV-lamp from above, when *E. coli* exists in a sample. On the other hand, blue color is observed in wells, when coliforms exist in a sample.

**[Determination]**

The total cell number is determined from the MPN calibration scale.

**[Purpose]**

A quantitative testing for *E. coli* (MPN testing using specific enzyme substrate-containing medium) was introduced as an official testing method for index bacterium of *Cryptosporidium* in tap water.

**[Warning and Direction for Use]**

**1. General precautions**

- 1) Read and follow precisely the warning and direction for use described on the package insert and/or label.
- 2) EC Blue solution is carefully poured into a slope of MPN plate to keep the solution from overflowing.
- 3) After incubation, MPN plate is gently treated to avoid a contamination among each of wells.
- 4) MPN plate is only for one-time use.

**2. Precautions for danger**

- 1) When media touched eyes or mouth, immediately wash with plenty of water, and consult a physician.
- 2) Manipulations with microorganisms always involve certain risks of laboratory-acquired infections. Manipulations should be done under the supervision of trained specialists with biohazard protection measures.

**3. Precautions for disposal of waste**

Any media, materials and apparatus after use must be sterilized by autoclaving or boiling water and then dispose as medical or industrial waste according to the Law on Waste Disposal and Cleaning.

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Note: The MPN calibration scale is based on the program developed in 1985 by Hurley & Roscoe. This program makes provision for the reliability coefficient to evaluate the reliability of data. In case that the combination of positive number indicates <9 of the reliability coefficient, MPN is valid. On the other hand, in case of >9, MPN may not be valid by some sort of error cause. MPN indicated in italic letter should be handled as a reference value because of >9.

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Shimadzu Diagnostics Corporation**