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** Easy Medium

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(9D17S)

Medium packaged into stick pouch

For detection of coliform bacteria

** Easy Medium

Desoxycholate Agar (Granule)

* Product is divided granulated media weighed and packaged into aluminum stick pouch for 400 mL. It can be used without weighing and dividing.

* Formula (in 45.0g for 1 liter)

Sodium Desoxycholate	1.0g
Peptone	10.0g
Ferric Ammonium Citrate	2.0g
Sodium Chloride	5.0g
Dipotassium Phosphate	2.0g
Lactose	10.0g
Neutral Red	0.033g
Agar	15.0g

pH 7.4 ± 0.1

* Directions

Open an aluminum pouch, put the medium (18.0g) into 400 mL of distilled water, mix well and heat to dissolve the medium at 100°C for 20 ~30minutes. Maintain the medium at about 50°C. Distribute about 15 mL into Petri dishes contained test samples in advance, and mix well to be solid. After being solid, pour about 5 mL to surface additionally as the second layer. Wait to be solid, and then incubate at 35 ± 1°C for 20 ± 2 hours. The medium does not require any sterilization.

Determination and Counting

Count red colonies grown on the plate as coliform bacteria number.

Purpose

Product is used for detection of coliform bacteria in foods including water, raw meats, fishes, raw vegetables and materials (vessel, chopping board and so on.) used together with the foods. Product is divided powder media weighed and packaged into aluminum stick pouch for 400 mL. It can be used without weighing and dividing.

Warning and Direction for Use

1. General precautions

- 1) Read and follow precisely the warning and direction for use described on the package insert and/or label.
- 2) Do not use the product after its expiry date. Quality of the product is not warranted after its shelf life.

- 3) Do not use the product with damaged aluminum pouch or that contains any foreign materials.
- 4) After cut aluminum pouch, recommend using it at once. Avoid reuse after storage of prepared media.
- 5) Add the medium to approximately one-third of the volume of distilled water. After thorough mixing, add the remainder of the water with care being taken to wash down the sided of the container.

2. Precautions for danger

- 1) When if medium touched eyes or mouth, immediately wash with plenty of water, and consult a physician.
- 2) Manipulations with microorganisms involve always certain risks of laboratory-acquired infections. Manipulations should be practiced under the supervision of key specialist with biohazard protection measures.

3. Precautions for disposal of waste

Any medium, materials and apparatus must be sterilized by autoclaving or boiling water after use, and then dispose as medical or industrial waste according to the Law on Waste Disposal and Cleaning. Also follow to local laws and regulations related to dispose such material after differentiation.

* Storage

Keep dry at room temperature (1 ~ 30°C).

Shelf Life

Three (3) years after manufacturing. Actual expiry date is printed on label of box and aluminum pouch.

** Package

Easy Medium Desoxycholate Agar

For 400 mL (18.0g) × 55 pouches Code 56204

** Further information

Customer Support Section,

Shimadzu Diagnostics Corporation

3-24-6, Ueno, Taito-ku, Tokyo 110-0005 Japan

TEL: +81-3-5846-5707

E-mail: contact@sd.c.shimadzu.co.jp

** Manufactured by

Shimadzu Diagnostics Corporation

3-24-6, Ueno, Taito-ku, Tokyo 110-0005 Japan

** Revised : April 2023

* Revised : April 2019
(9D17S)