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(0257S)

Medium packaged into stick pouch
For isolation of *Staphylococcus aureus*

** Easy Medium

Salt Egg Yolk Agar Base (Granule)

Product is divided powder media weighed and packaged into aluminum stick pouch for 200 mL. It can be used without weighing and dividing.

Formula (in 115.0g for 1 liter)

Heart Extract Powder	5.0g
Peptone	10.0g
Sodium Chloride	70.0g
Lithium Chloride	5.0g
Mannitol	10.0g
Agar	15.0g

pH 7.3 ± 0.2

Directions

Put one pouch (23.0 g) into 200 mL of distilled water, mix well and heat to dissolve the medium. Sterilize by autoclaving at 121°C for 15 minutes. Cool down to approximately 50°C. Add an egg yolk suspension in a ratio of 10% with thorough mixing. Distribute about 20 mL into Petri dishes. Dry the surface of the plate before use. Inoculate 0.1 mL of test sample, and incubate it at 35±1°C for 48 hours.

[Preparation of Egg Yolk Solution]

Fresh egg is sterilized by alcohol, put the egg yolk alone to sterilized cylinder aseptically, add equivalent volume of 10% sterilized sodium chloride solution, and mix well to have homogeneous solution.

Determination and Counting

Staphylococcus aureus generate yellow/lemon (sometimes white) circle colony with pearl turbid ring (halo) around the colony by egg yolk reaction.

Purpose

Product is used for isolation of *Staphylococcus aureus* in foods or from environment of food production and cooking (apparatus, devices or workers fingers).

Staphylococcus aureus can be distinguished by egg yolk reaction.

Warning and Direction for Use

1. General precautions

- 1) Read and follow precisely the warning and direction for use described on the package insert and/or label.
- 2) Do not use the product after its expiry date. Quality of the product is not warranted after its shelf life.
- 3) Do not use the product with damaged aluminum pouch or that contains any foreign materials.
- 4) After cut aluminum pouch, recommend using it at once. Avoid reuse after storage of prepared media.

2. Precautions for danger

- 1) When if medium touched eyes or mouth, immediately wash with plenty of water, and consult a physician.
- 2) Manipulations with microorganisms involve always certain risks of laboratory-acquired infections. Manipulations should be practiced under the supervision of key specialist with biohazard protection measures.

3. Precautions for disposal of waste

Any medium, materials and apparatus must be sterilized by autoclaving or boiling water after use, and then dispose as medical or industrial waste according to the Law on Waste Disposal and Cleaning. Also follow to local laws and regulations related to dispose such material after differentiation.

Storage

Keep dry at room temperature.

Shelf Life

Three (3) years after manufacturing. Actual expiry date is printed on label of box and aluminum pouch.

** Package

Easy Medium Salt Egg Yolk Agar Base (Granule)

For 200 mL (23.0g) x 44 pouches Code 56216

*** Further information

Customer Support Section,

Shimadzu Diagnostics Corporation

3-24-6, Ueno, Taito-ku, Tokyo 110-0005 Japan

TEL: +81-3-5846-5707 / FAX: +81-3-5846-5629

E-mail: contact@cdc.shimadzu.co.jp

** * Manufactured by

Shimadzu Diagnostics Corporation

3-24-6, Ueno, Taito-ku, Tokyo 110-0005 Japan

** Revised : April 2023

* Revised : January 2018
(0257S)