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(0A06S)

Medium packaged into stick pouch  
For detection of *E. coli* and coliform bacteria by colorimetric enzyme substrate

## \*\* Easy Medium XM-G Agar (Granule)

\* Product is divided granulated media weighed and packaged into aluminum stick pouch for 400 mL. It can be used without weighing and dividing.

### Formula (in 39.3g for 1 liter)

|   |       |
|---|-------|
| Peptone   | 10.0g |
| Sodium Pyruvate   | 1.0g  |
| L-Tryptophan  | 1.0g  |
| D-Sorbitol  | 1.0g  |
| Sodium Chloride   | 5.0g  |
| Monopotassium Phosphate   | 2.2g  |
| Dipotassium Phosphate   | 2.7g  |
| Potassium Nitrate   | 1.0g  |
| Sodium Lauryl Sulfate   | 0.2g  |
| 5-Bromo-4-Chloro-3-Indolyl- $\beta$ -<br>D-Glucuronide(X-GLUC)            | 0.1g  |
| 5-Bromo-6-Chloro-3-Indolyl- $\beta$ -<br>D-Galactopyranoside(MAGENTA-GAL) | 0.1g  |
| Agar  | 15.0g |
| pH 7.0 $\pm$ 0.2  |       |

### \* Directions

Put one pouch (15.7 g) into 400 mL of distilled water, mix well and sterilize by autoclaving at 121°C for 15 minutes. Cool down to around 50°C, and distribute it into Petri dish for use. For liquid specimen, the sterilized medium is distributed into Petri dish contained specimen dispensed in advance to make pour culture method. Incubate it at 35  $\pm$  1°C for 20  $\pm$  2 hours, and then make judgment.

### Determination and Counting

After incubation, blue to bluish purple colonies are judged as *E. coli*, and pink to red colonies are judged as coliform bacteria. *E. coli* has specific enzyme of  $\beta$ -Glucuronidase that reacts with substrate of  $\beta$ -GLUC to develop blue color. Coliform bacteria has specific enzyme of  $\beta$ -Galactosidase that reacts with substrate of MAGENTA-GAL to develop red color. Since *E. coli* has both enzymes, and then mixed color of bluish purple to purple is also developed. O-157 does not have the enzyme of  $\beta$ -Glucuronidase, and then is detected as coliform bacteria. For specimen containing *Lactobacillus*, red color may be developed in some cases by its enzyme. On such cases, confirmation test like gram staining is recommended to make in order to distinguish.

### Purpose

Product is used for detection of *E. coli* and coliform bacteria in foods, refreshing drinks, and wiped materials.

### Warning and Direction for Use

#### 1. General precautions

- 1) Read and follow precisely the warning and direction for use described on the package insert and/or label.
- 2) Do not use the product after its expiry date. Quality of the product is not warranted after its shelf life.
- 3) Do not use the product with damaged aluminum pouch or that contains any foreign materials.
- 4) After cut aluminum pouch, recommend using it at once. Avoid reuse after storage of prepared media.

#### 2. Precautions for danger

- 1) When if medium touched eyes or mouth, immediately wash with plenty of water, and consult a physician.
- 2) Manipulations with microorganisms involve always certain risks of laboratory-acquired infections. Manipulations should be practiced under the supervision of key specialist with biohazard protection measures.

#### 3. Precautions for disposal of waste

Any medium, materials and apparatus must be sterilized by autoclaving or boiling water after use, and then dispose as medical or industrial waste according to the Law on Waste Disposal and Cleaning. Also follow to local laws and regulations related to dispose such material after differentiation.

### Storage

Keep dry at room temperature.

### Shelf Life

Three (3) years after manufacturing. Actual expiry date is printed on label of box and aluminum pouch.

### Package

Easy Medium XM-G Agar

For 400 mL (15.7g) x 30 pouches ..... Code 56228

### \*\* Further information

#### Customer Support Section,

#### Shimadzu Diagnostics Corporation

3-24-6, Ueno, Taito-ku, Tokyo 110-0005 Japan

TEL: +81-3-5846-5707

E-mail: [contact@sdcs.shimadzu.co.jp](mailto:contact@sdcs.shimadzu.co.jp)

\*\* Manufactured by

**Shimadzu Diagnostics Corporation**

3-24-6, Ueno, Taito-ku, Tokyo 110-0005 Japan

\*\* Revised : April 2023

\* Revised : January 2020  
(0A06S)