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(9H03S)

Medium packaged into stick pouch
Medium with chromogenic substrate for detection of *Staphylococcus aureus*

** Easy Medium

X-SA Agar (Granule)

Product is divided granulated powder media weighed and packaged into aluminum stick pouch for 200 mL. It can be used without weighing and dividing.

Formula (in 51.0g for 1 liter)

Peptone	13.0g
Beef Extract	3.0g
Lithium Sulfate	10.0g
Mannitol	10.0g
Agar	14.0g
Selective agent	0.7g
Chromogenic Substrate	0.3g
pH 7.3 ± 0.2	

*Directions

Put one pouch (10.2 g) into 200 mL of distilled water, mix well and heat to dissolve the medium at 100°C for 20-30 minutes. Maintain the medium at about 50°C. Distribute about 20 mL of medium into Petri dishes for use. (Do not sterilize by autoclaving.)

Inoculate specimen, incubate the medium at 35 - 37°C for 22 - 24 hours, and after incubation make judgment. When incubating longer than 24 hours, other *Staphylococcus* than *Staphylococcus aureus* may also generate light blue colony. Then make judgment within 24 hours in order to avoid misjudgment.

Determination and Counting

After incubation, blue and light blue colonies are counted as *Staphylococcus aureus*.

Coagulase Negative *Staphylococcus* (CNS) like *Staphylococcus epidermidis* generate white or pale light blue minimal colony, and so they are easily distinguished. Other tests like gram staining or coagulase test are recommended to make in order to confirm *Staphylococcus aureus*.

Although other bacteria than *Staphylococcus* are inhibited to grow by selective agent, some genus *Bacillus* may grow and generate light blue colony. They generate flat and lusterless colony, and then distinguish by gram staining.

Purpose

Product is used for detection of *Staphylococcus aureus* in foods or wiped materials.

Warning and Direction for Use

1. General precautions

- 1) Read and follow precisely the warning and direction for use described on the package insert and/or label.
- 2) Do not use the product after its expiry date. Quality of the product is not warranted after its shelf life.
- 3) Do not use the product with damaged aluminum pouch or that contains any foreign materials.
- 4) After cut aluminum pouch, recommend using it at once. Avoid reuse after storage of prepared media.

2. Precautions for danger

- 1) When if medium touched eyes or mouth, immediately wash with plenty of water, and consult a physician.
- 2) Manipulations with microorganisms involve always certain risks of laboratory-acquired infections. Manipulations should be practiced under the supervision of key specialist with biohazard protection measures.

3. Precautions for disposal of waste

Any medium, materials and apparatus must be sterilized by autoclaving or boiling water after use, and then dispose as medical or industrial waste according to the Law on Waste Disposal and Cleaning. Also follow to local laws and regulations related to dispose such material after differentiation.

Storage

Keep dry at room temperature(1~30°C).

Shelf Life

Two (2) years after manufacturing. Actual expiry date is printed on label of box and aluminum pouch.

** Package

Easy Medium X-SA Agar (Granule)

For 200 mL (10.2g) x 30 pouches Code 56230

** * Further information

Customer Support Section,

Shimadzu Diagnostics Corporation

3-24-6, Ueno, Taito-ku, Tokyo 110-0005 Japan

TEL: +81-3-5846-5707 / FAX: +81-3-5846-5629

E-mail: contact@sd.c.shimadzu.co.jp

** * Manufactured by

Shimadzu Diagnostics Corporation

3-24-6, Ueno, Taito-ku, Tokyo 110-0005 Japan

** Revised : April 2023

* Revised : September 2019
(9H04S)